

Sep Sci

HHC
Helvetica Health Care

BIOFIRE SPOTFIRE CONTROLS

NEW!

**NATRSTM-BIO - NATtrol Respiratory/
Sore Throat (R/ST) Verification Panel
Mini 7 x 0.6mL; 4 x 1.8 mL in stock**

NATRSTM-BIO

NATtrol Respiratory/Sore Throat (R/ST) Verification Panel Mini (NATRSTM-BIO) can be used to develop procedures for verification of performance of molecular assays that detect the presence of bacterial and viral nucleic acids. The panel can also be used for training and evaluation of operator proficiency.

NATtrol Respiratory/Sore Throat Verification Panel Mini contains purified, intact bacteria and viruses that have been chemically modified to render them non-infectious and refrigerator stable. The panel members are supplied in a proprietary matrix.



Quality Control Panels



Validation-Verification

Contact **Sep Sci** for more information.

Separation Scientific

Tel: 011 794 4117 | Email: Info@sepsci.co.za | Web: www.sepsci.co.za