



Encephalitis Icon

Now offering NEW heat inactivated culture fluids for tick-borne Encephalitis virus



Encephalitis is inflammation of the active tissues of the brain caused by an infection or an autoimmune response.



Tick-Borne Encephalitis Virus (Strain: Neudorfl) and (Strain: Vaccine) is an enveloped capsid virus 40-60 nm in diameter and is a positive-sense single-stranded RNA virus. Each frozen aliquot contains 1 mL of heat inactivated viral culture fluid. The pre-inactivation titer was determined from an infectious aliquot. Viral inactivation is verified after heat inactivation by the absence of viral growth in tissue culture-based infectivity assays.

Heat inactivated viral culture fluids are sold as consumable testing materials. The suitability and performance characteristics should be determined by your laboratory for each intended usage.

Contact **Sep Sci** for more information.



Separation Scientific

Tel: 011 794 4117 | Email: Info@sepsci.co.za | Web: www.sepsci.co.za